

## Product datasheet for **RC202168L3V**

### Calpastatin (CAST) (NM\_001042444) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Calpastatin (CAST) (NM_001042444) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Calpastatin
Synonyms:	BS-17; PLACK
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001042444
ORF Size:	2001 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202168).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001042444.1</a>
RefSeq Size:	4366 bp
RefSeq ORF:	2004 bp
Locus ID:	831
UniProt ID:	<a href="#">P20810</a>
Cytogenetics:	5q15
MW:	71.9 kDa


[View online »](#)

**Gene Summary:**

The protein encoded by this gene is an endogenous calpain (calcium-dependent cysteine protease) inhibitor. It consists of an N-terminal domain L and four repetitive calpain-inhibition domains (domains 1-4), and it is involved in the proteolysis of amyloid precursor protein. The calpain/calpastatin system is involved in numerous membrane fusion events, such as neural vesicle exocytosis and platelet and red-cell aggregation. The encoded protein is also thought to affect the expression levels of genes encoding structural or regulatory proteins. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jun 2010]