

Product datasheet for **RC202181L4V**

BAHD1 (NM_014952) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	BAHD1 (NM_014952) Human Tagged ORF Clone Lentiviral Particle
Symbol:	BAHD1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_014952
ORF Size:	2340 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202181).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_014952.2 , NP_055767.2
RefSeq Size:	4605 bp
RefSeq ORF:	2343 bp
Locus ID:	22893
UniProt ID:	Q8TBE0
Cytogenetics:	15q15.1
MW:	84.7 kDa



[View online »](#)

Gene Summary:

Heterochromatin protein that acts as a transcription repressor and has the ability to promote the formation of large heterochromatic domains. May act by recruiting heterochromatin proteins such as CBX5 (HP1 alpha), HDAC5 and MBD1. Represses IGF2 expression by binding to its CpG-rich P3 promoter and recruiting heterochromatin proteins. At specific stages of Listeria infection, in complex with TRIM28, corepresses interferon-stimulated genes, including IFNL1, IFNL2 and IFNL3.[UniProtKB/Swiss-Prot Function]