

Product datasheet for **RC202185A1V**

Human hHR23b (RAD23B) (NM_002874) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human hHR23b (RAD23B) (NM_002874) AAV Particle
Tag:	Myc-DDK
Symbol:	hHR23b
Synonyms:	HHR23B; HR23B; P58
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC202185 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCAGGTCACCCTGAAGACCCTCCAGCAGCAGACCTTCAAGATAGACATTGACCCCGAGGAGACGGTGA
AAGCACTGAAAGAGAAGATTGAATCTGAAAAGGGGAAAGATGCCTTCCAGTAGCAGGTCAAAAATTAAT
TTATGCAGGCAAAATCCTCAATGATGATACTGCTCTCAAAGAATATAAAATTGATGAGAAAACTTTGTG
GTGGTTATGGTGACCAAAACCAAGCAGTGTCCACACCAGCAGCTACAACCTCAGCAGTCAGCTCCTG
CCAGCACTACAGCAGTTACTTCTCCACCACCACAAGTGTGGCTCAGGCTCCAACCCCTGTCCCTGCCTT
GGCCCCCACTTCCACACCTGCATCCATCACTCCAGCATCAGCGACAGCATCTTCTGAACCTGCACCTGCT
AGTGCAGCTAAACAAGAGAAGCCTGCAGAAAAGCCAGCAGAGACACCAAGTGGCTACTAGCCCAACAGCAA
CTGACAGTACATCGGGTGATTCTTCTCGGTCAAACCTTTTGAAGATGCAACGAGTGCCTTGTGACGGG
TCAGTCTTACGAGAATATGGTAACTGAGATCATGTCAATGGGCTATGAACGAGAGCAAGTAATTGCAGCC
CTGAGAGCCAGTTTCAACAACCTGACAGAGCAGTGGAGTATCTTTAATGGGAATCCCTGGAGATAGAG
AAAGTCAGGCTGTGGTTGACCCCTCAAGCAGCTAGTACTGGGTTCTCAGTCTTCAGCAGTGGCTGC
AGCTGCAGCAACTACGACAGCAACAACATAACAACAAGTCTGGAGGACATCCCCTTGAATTTTACGG
AATCAGCCTCAGTTTCAACAGATGAGACAAATTATTCAGCAGAATCCTTCTTGTCTCCAGCGTTACTAC
AGCAGATAGGTGAGAGAATCCTCAATTACTTCAGCAAATTAGCCAACACAGGAGCATTATTCAGAT
GTTAAATGAACAGTTCAAGAAGCTGGTGGTCAAGGAGGAGGAGTGGAGGTGGCAGTGGAGGAATTGCA
GAAGCTGGAAGTGGTCATATGAACTACATTCAAGTAACACCTCAGGAAAAAGAAGCTATAGAAAGGTTAA
AGGCATTAGGATTTCTGAAGGACTTGTGATACAAGCGTATTTTGCTTGTGAGAAGAATGAGAATTTGGC
TGCCAATTTTCTTCTACAGCAGAAGCTTTGATGAAGAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC202185 protein sequence Red=Cloning site Green=Tags(s) MQVTLKTLQQQTFKIDIDPEETVKALKEKIESEKGDAPFPVAGQKLIYAGKILNDDTALKEYKIDEKNFV VVMVTKPKAVSTPAPATTQQSAPASTTAVTSSTTTTVAQPTPVPALAPTSTPASITPASATASSEPAPA SAAKQEKPAEKPAETPVATSPTATDSTSGDSSRSNLFEDATSALVTGQSYENMVTEIMSMGYEREQVIAA LRASFNNPDRAVEYLLMGIPGDRESQAVDPPQAASTGVPQSSAVAAAAATTTATTTTSSGGHPLEFLR NQPFQQMRQIIQQNPSSLPALLLQIGRENPLLQQISQHQEHFIQMLNEPVQEAGGQGGGGGGGGGIA EAGSGHMNYIQVTPQEKEAIERLKALGFPEGLVIQAYFACEKNENLAANFLLQQNFDED TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_002874
ORF Size:	1227 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_002874.3
RefSeq Size:	4167 bp
RefSeq ORF:	1230 bp
Locus ID:	5887
UniProt ID:	P54727
Cytogenetics:	9q31.2
MW:	43.2 kDa