

Product datasheet for RC202188A1V

Human MAGP1 (MFAP2) (NM_017459) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human MAGP1 (MFAP2) (NM_017459) AAV Particle
Tag: Myc-DDK
Symbol: MAGP1
Synonyms: MAGP; MAGP-1; MAGP1
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC202188 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC

ATGAGAGCTGCCTACCTCTTCTGCTATTCTGCCTGCAGGCTTGCTGGCTCAGGGCCAGTATGACCTGG
 ACCCGCTGCCGCCGTTCCCTGACCACGTCCAGTACACCCACTATAGCGACCAGATCGACAACCCAGACTA
 CTATGATTATCAAGAGGTGACTCCTCGGCCCTCCGAGGAACAGTTCAGTTCAGTCCAGCAGCAAGTC
 CAACAGGAAGTCATCCAGCCCCAACCCAGAACAGGAAATGCAGAGCTGGAGCCACAGAGCCTGGGC
 CTCTTGACTGCCGTGAGGAACAGTACCCGTGCACCCGCCCTCTACTCCATACACAGGCCTTGCAACAGTG
 TCTCAACGAGGTCTGCTTCTACAGCCTCCGCCGTGTGTACGTCAATTAACAAGGAGATCTGTGTTCTGTACA
 GTGTGTGCCACGAGGAGCTCCTCCGAGCTGACCTCTGTGCGGACAAGTTCTCCAAATGTGGCGTGATGG
 CCAGCAGCGGCCTGTGCAATCCGTGGCGGCCTCTGTGCCAGGAGCTGTGGGAGCTGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC202188 protein sequence
 Red=Cloning site Green=Tags(s)

MRAAYLFLLFLPAGLLAQGQYDLDPFPDPHVQYTHYSQIDNPDIYDYQEVTPRPSEEQFQFQSQQQV
 QQEVIPAPTPEPGNAELEPTEPGPLDCREEQYPCTRLYSIHRPCKQCLNEVCFYSLRRVYVINKEICVRT
 VCAHEELLRADLCRDKFSKCGVMASSGLCQSVAAASCARSCGSC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2


[View online »](#)

ACCN:	NM_017459
ORF Size:	549 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_017459.1</u>
RefSeq Size:	1105 bp
RefSeq ORF:	552 bp
Locus ID:	4237
UniProt ID:	<u>P55001</u>
Cytogenetics:	1p36.13
MW:	20.8 kDa