

Product datasheet for **RC202197L4V**

RAB24 (NM_130781) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RAB24 (NM_130781) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RAB24
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_130781
ORF Size:	609 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202197).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_130781.2 , NP_570137.2
RefSeq Size:	1546 bp
RefSeq ORF:	612 bp
Locus ID:	53917
UniProt ID:	Q969Q5
Cytogenetics:	5q35.3
Domains:	ras, RAS, RHO, RAB
Protein Families:	Druggable Genome
MW:	23.1 kDa



[View online »](#)

Gene Summary:

RAB24 is a small GTPase of the Rab subfamily of Ras-related proteins that regulate intracellular protein trafficking (Oikkonen et al., 1993 [PubMed 8126105]).[supplied by OMIM, Aug 2009]