

## Product datasheet for **RC202200L2V**

### **FKBP25 (FKBP3) (NM\_002013) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	FKBP25 (FKBP3) (NM_002013) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FKBP25
Synonyms:	FKBP-3; FKBP-25; FKBP25; PPIase
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_002013
ORF Size:	672 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202200).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_002013.2</a>
RefSeq Size:	1353 bp
RefSeq ORF:	675 bp
Locus ID:	2287
UniProt ID:	<a href="#">Q00688</a>
Cytogenetics:	14q21.2
Domains:	FKBP
Protein Families:	Druggable Genome, Stem cell - Pluripotency



[View online »](#)

**MW:** 25.2 kDa

**Gene Summary:** The protein encoded by this gene is a member of the immunophilin protein family, which play a role in immunoregulation and basic cellular processes involving protein folding and trafficking. This encoded protein is a cis-trans prolyl isomerase that binds the immunosuppressants FK506 and rapamycin, as well as histone deacetylases, the transcription factor YY1, casein kinase II, and nucleolin. It has a higher affinity for rapamycin than for FK506 and thus may be an important target molecule for immunosuppression by rapamycin. [provided by RefSeq, Sep 2008]