

Product datasheet for **RC202206A1V**

Human GALT (NM_000155) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human GALT (NM_000155) AAV Particle
Tag:	Myc-DDK
Symbol:	GALT
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC202206 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTCGCGCAGTGAACCGATCCTCAGCAACGCCAGCAGGCGTCAGAGGCGGACGCCGAGCAGCAACCT
TCCGGGCAAACGACCATCAGCATATCCGCTACAACCCGCTGCAGGATGAGTGGGTGCTGGTGTCAGCTCA
CCGCATGAAGCGGCCCTGGCAGGTCAGTGGAGCCCCAGCTTCTGAAGACAGTGCCCCGCCATGACCT
CTCAACCTCTGTGTCTGGGGCCATCCGAGCCAACGGAGAGGTGAATCCCCAGTACGATAGCACCTTCC
TGTTTGACAACGACTTCCAGCTCTGCAGCCTGATGCCCCAGTCCAGGACCCAGTGATCATCCCTTTT
CCAAGCAAAGTCTGCTCGAGGAGTCTGTAAGGTCATGTGCTTCCACCCCTGGTCGGATGTAACGCTGCCA
CTCATGTCCGTCCTGAGATCCGGGCTGTTGTTGATGCATGGGCTCAGTCACAGAGGAGCTGGGTGCC
AGTACCCTTGGGTGCAGATCTTTGAAAACAAAGGTGCCATGATGGGCTGTTCTAACCCCCACCCCACTG
CCAGGTATGGGCCAGCAGTTTCTGCCAGATATTGCCAGCGTGAGGAGCGATCTCAGCAGGCCTATAAG
AGTCAGCATGGAGAGCCCTGCTAATGGAGTACAGCCGCCAGGAGCTACTCAGGAAGGAACGTCTGGTCC
TAACAGTGAGCACTGGTTAGTACTGGTCCCCTTCTGGCAACATGGCCCTACCAGACACTGCTGCTGCC
CCGTCGGCATGTGCGGCGGCTACCTGAGCTGACCCCTGCTGAGCGTGATGATCTAGCCTCCATCATGAAG
AAGCTCTTGACCAAGTATGACAACCTCTTTGAGACGTCCTTTCCCTACTCCATGGGCTGGCATGGGGCTC
CCACAGGATCAGAGGCTGGGGCCAACCTGGGACCATTTGGCAGCTGCACGCTCATTACTACCTCCGCTCCT
GCGCTCTGCCACTGTCCGAAATTCATGGTTGGCTACGAAATGCTTGCTCAGGCTCAGAGGGACCTCACC
CCTGAGCAGGCTGCAGAGAGACTAAGGGCACTTCTGAGGTTTCATTACCACCTGGGGCAGAAGGACAGGG
AGACAGCAACCATCGCC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC202206 protein sequence Red=Cloning site Green=Tags(s) MSRSGTDPQQRQQASEADAAAATFRANDHQHIRYNPLQDEWVLVSAHRMKRPWQGQVEPQLLKTVPRHDP LNPLCPGAIRANGEVNPQYDSTFLFDNDFPALQPDAPSPGSDHPLFQAKSARGVCKVMCFHPWSDVTLF LMSVPEIRAVVDASVTEELGAQYPWVQIFENKGAMMGCSNPHPHCQVWASSFLPDIAQREERSQQAYK SQHGEPDLLMEYSRQELLRKERLVL TSEHWLVLVPFWATWPYQTLTLLPRRHVRRLELTPAERDDLASIMK KLLTKYDNLFETSFYPSMGWHGAPTGSEAGANWDHWQLHAHYYPPLRSATVRKFMVGYEMLAQQRDLT PEQAAERLRALPEVHYHLGQKDRETATIA TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_000155
ORF Size:	1137 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_000155.2</u>
RefSeq Size:	1407 bp
RefSeq ORF:	1140 bp
Locus ID:	2592
UniProt ID:	<u>P07902</u>
Cytogenetics:	9p13.3
MW:	43.4 kDa