

Product datasheet for **RC202215**

ZMAT5 (NM_019103) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: ZMAT5 (NM_019103) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: ZMAT5
Synonyms: SNRNP20; U11/U12-20K; ZC3H19
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC202215 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGGAAGCGATACTTCTGTGACTACTGCGACCCTCTTCCAGGACAACCTCCACAACCGCAAGAAGC
ACCTGAACGGGCTGCAGCACCTCAAGGCCAAGAAGGTCTGGTACGACATGTTCCGAGATGCAGCTGCCAT
CTTGCTGGATGAGCAGAACAAGCGGCCCTGCAGGAAGTTCTACTGACAGGCCAGTGCAGCTTTGGCTCC
AACTGCAGATTTCCACATGTCAGAGCGAGACCTGCAGGAGCTGAGCATCCAGGTGGAGGAGGAGAGGC
GAGCCAGGGAGTGGCTACTAGATGCTCCTGAGCTCCCCGAGGGCCATCTGGAGGACTGGCTGGAGAAGAG
AGCCAAGCGGCTGAGCTCAGCCCAAGTAGCAGGGCTGAACCCATCAGAACCACTGTCTTCCAGTACCCC
GTGGGCTGGCCACCAGTTCAGGAGCTGCCTCCATCCCTGCGGGCACCCCCACCTGGGGGTGGCCTCTGC
AGCCAGAGTCCAGTGGGGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC202215 protein sequence
Red=Cloning site Green=Tags(s)

MGKRYFCDYCDRSFQDNLHNRKKHLNGLQHLKAKKVWYDMFRDAAAILLDEQNKRPCRKFLLTGQCDFGS
NCRFSHMSERDLQELSIQVEEERRAREWLLDAPELPEGHLEDWLEKRAKRLSSAPSSRAEPIRTTVFYQP
VGWPPVQELPPSLRAPPFGGWPLQPRVQWG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Chromatograms: https://cdn.origene.com/chromatograms/mk6308_d04.zip



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_019103

ORF Size: 510 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_019103.2](#), [NP_061976.1](#)

RefSeq Size: 1051 bp

RefSeq ORF: 513 bp

Locus ID: 55954


UniProt ID: [Q9UDW3](#)

Cytogenetics: 22q12.2


Domains: zf-CCCH

MW: 20 kDa


Product images:



Circular map for RC202215



Western blot validation of overexpression lysate (Cat# [LY423980]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with [RC218362] using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified ZMAT5 protein (Cat# [TP302215]). The protein was produced from HEK293T cells transfected with ZMAT5 cDNA clone (Cat# RC202215) using MegaTran 2.0 (Cat# [TT210002]).