

## Product datasheet for **RC202216L2V**

### EMP3 (NM\_001425) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	EMP3 (NM_001425) Human Tagged ORF Clone Lentiviral Particle
Symbol:	EMP3
Synonyms:	YMP
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_001425
ORF Size:	489 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202216).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001425.1</a>
RefSeq Size:	850 bp
RefSeq ORF:	492 bp
Locus ID:	2014
UniProt ID:	<a href="#">P54852</a>
Cytogenetics:	19q13.33
Domains:	PMP22_Claudin
Protein Families:	Transmembrane



[View online »](#)

**MW:** 18.4 kDa

**Gene Summary:** The protein encoded by this gene belongs to the PMP-22/EMP/MP20 family of proteins. The protein contains four transmembrane domains and two N-linked glycosylation sites. It is thought to be involved in cell proliferation, cell-cell interactions and function as a tumor suppressor. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2015]