

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# Product datasheet for RC202253L4V

### MYD88 (NM\_002468) Human Tagged ORF Clone Lentiviral Particle

#### **Product data:**


Product Type:	Lentiviral Particles
Product Name:	MYD88 (NM_002468) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MYD88
Synonyms:	IMD68; MYD88D
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_002468
ORF Size:	888 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202253).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<u>NM 002468.2</u>
RefSeq Size:	2862 bp
RefSeq ORF:	891 bp
Locus ID:	4615
UniProt ID:	<u>Q99836</u>
Cytogenetics:	3p22.2
Domains:	TIR, DEATH
Protein Families:	Druggable Genome



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	MYD88 (NM_002468) Human Tagged ORF Clone Lentiviral Particle – RC202253L4V
Protein Pathways	Apoptosis, Toll-like receptor signaling pathway
MW:	33.2 kDa
Gene Summary:	This gene encodes a cytosolic adapter protein that plays a central role in the innate and adaptive immune response. This protein functions as an essential signal transducer in the interleukin-1 and Toll-like receptor signaling pathways. These pathways regulate that activation of numerous proinflammatory genes. The encoded protein consists of an N- terminal death domain and a C-terminal Toll-interleukin1 receptor domain. Patients with defects in this gene have an increased susceptibility to pyogenic bacterial infections. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Feb 2010]

## Product images:



[RC202253L4] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC202253L4V particle to overexpress human MYD88-mGFP fusion protein.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US