

## Product datasheet for **RC202273L2V**

### **ARHGEF19 (NM\_153213) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	ARHGEF19 (NM_153213) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ARHGEF19
Synonyms:	WGEF
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_153213
ORF Size:	2406 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202273).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_153213.2</a>
RefSeq Size:	3044 bp
RefSeq ORF:	2409 bp
Locus ID:	128272
UniProt ID:	<a href="#">Q8IW93</a>
Cytogenetics:	1p36.13
MW:	89.2 kDa



[View online »](#)

**Gene Summary:**

Guanine nucleotide exchange factors (GEFs) such as ARHGEF19 accelerate the GTPase activity of Rho GTPases (see RHOA, MIM 165390).[supplied by OMIM, Dec 2008]