

## Product datasheet for **RC202285L4V**

### Apoptosis enhancing nuclease (AEN) (NM\_022767) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Apoptosis enhancing nuclease (AEN) (NM_022767) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Apoptosis enhancing nuclease
Synonyms:	ISG20L1; pp12744
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_022767
ORF Size:	975 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202285).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_022767.2</a> , <a href="#">NP_073604.2</a>
RefSeq Size:	3134 bp
RefSeq ORF:	978 bp
Locus ID:	64782
UniProt ID:	<a href="#">Q8WTP8</a>
Cytogenetics:	15q26.1
Domains:	EXOIII



[View online »](#)

**MW:** 36.4 kDa

**Gene Summary:** Exonuclease with activity against single- and double-stranded DNA and RNA. Mediates p53-induced apoptosis. When induced by p53 following DNA damage, digests double-stranded DNA to form single-stranded DNA and amplifies DNA damage signals, leading to enhancement of apoptosis.[UniProtKB/Swiss-Prot Function]