

## Product datasheet for RC202297A1V

### Human LRTOMT (NM\_145309) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Human LRTOMT (NM\_145309) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** LRTOMT  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >RC202297 ORF sequence  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCGCGATCGCC

ATGAACAAACGGGACTATATGAACACTTCGGTACAGGAGCCCCCTCTTGACTACTCCTTCAGAAGCATCC  
 ACGTCATTCAAGATCTGGTAAATGAGGAGCCAAGGACAGGACTACGACCACTGAAGCGTTCAAAGTCGGG  
 GAAATCACTGACCCAGTCCCTGTGGCTGAATAACAATGTTCTCAATGATCTGAGAGACTTCAACCAGGTG  
 GCTTCACAGCTGTTGGAGCACCCAGAGAACCTGGCCTGGATCGACCTGTCCTTTAATGACCTGACTTCCA  
 TTGACCCTGTCCTAACAACTTTCTTCAACCTGAGTGTCTCTATCTTCACGGCAACAGCATCCAGCGCCT  
 GGGGGAGGTGAATAAGCTGGCTGCTTCCCTCGGCTCCGTAGCCTGACACTCCATGGGAACCCCATGGAG  
 GAAGAGAAAGGGTATAGGCAATATGTGCTGTGCACCCTGTCCCGTATCACCACGTTGACTTCAGTGGGG  
 TCACCAAAGCAGACCGCACACAGCTGAAGTCTGGAAACGCATGAACATCAAGCCCAAGAAGGCCTGGAC  
 CAAGCAGAATACACTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC202297 protein sequence  
 Red=Cloning site Green=Tags(s)

MNKRDMNTSVQEPPLDYSFRSIHVIQDLVNEEPRTGLRPLKRSKSGKSLTQSLWLNNNVLNDRDFNQV  
 ASQLEHPENLAWIDLTFNDLTSIDPVLTTFFNLVLYLHGNSIQRLGEVNKLAVLPRLRSLTLHGPNME  
 EEKGYRQYVLCITLSRITTFDFSGVTADRTTAEVWKRNNIKPKKAWTKQNTL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Human  
**Serotype:** AAV-2


[View online »](#)

|                      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>ACCN:</b>         | NM_145309                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>ORF Size:</b>     | 576 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Buffer:</b>       | PBS with 0.001% Pluronic F68                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Stability:</b>    | AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus. |
| <b>RefSeq:</b>       | <u><a href="#">NM_145309.1</a></u>                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>RefSeq Size:</b>  | 2658 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>RefSeq ORF:</b>   | 579 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Locus ID:</b>     | 220074                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>UniProt ID:</b>   | <u><a href="#">Q96E66</a></u>                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Cytogenetics:</b> | 11q13.4                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>MW:</b>           | 22.2 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |