

Product datasheet for **RC202305L4V**

MED30 (NM_080651) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	MED30 (NM_080651) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MED30
Synonyms:	MED30S; THRAP6; TRAP25
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_080651
ORF Size:	534 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202305).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_080651.1
RefSeq Size:	994 bp
RefSeq ORF:	537 bp
Locus ID:	90390
UniProt ID:	Q96HR3
Cytogenetics:	8q24.11
Protein Families:	Druggable Genome, Transcription Factors
MW:	20.3 kDa



[View online »](#)

Gene Summary:

The multiprotein TRAP/Mediator complex facilitates gene expression through a wide variety of transcriptional activators. MED30 is a component of this complex that appears to be metazoan specific (Baek et al., 2002 [PubMed 11909976]).[supplied by OMIM, Nov 2010]