

Product datasheet for **RC202326A1V**

Human DIMT1L (DIMT1) (NM_014473) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human DIMT1L (DIMT1) (NM_014473) AAV Particle
Tag:	Myc-DDK
Symbol:	DIMT1L
Synonyms:	DIM1; DIMT1L; HSA9761; HUSSY5
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC202326 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCCGAAGGTCAAGTCGGGGGCCATCGGCCGCCGCCGGCGGCAGGAGCAGCGCCGGGAGCTGAAGA
GCGCTGGAGGACTCATGTTCAACACGGGGATTGGGCAGCACATTTGAAAAATCCTCTCATTATTAACAG
CATTATCGATAAGGCTGCCTTAAGACCAACTGATGTAGTGCTGGAAGTTGGACCTGGAAGTGGCAACATG
ACTGTAAAGTTGTTAGAAAAGGCAAAAAGTTGTTGCTTGTAAGTTGACCAAGGCTAGTAGCTGAAC
TTCACAAAAGAGTTCAAGGCACGCTGTGGCCAGCAAACTTCAAGTACTGGTGGGTGATGTGCTGAAAAC
AGATTTGCCATTCTTTGATACTTGTGTGGCAAATTTGCCCTTATCAGATCTCTTCGCCTTTTGTCTTCAAG
CTGTTGCTACATCGACCTTTTTTCAGGTGTGCTATACTTATGTTTCAAAGAGAATTTGCCCTCCGACTGG
TTGCAAAACCTGGAGATAAGTTATACTGCAGACTCTCAATTAATACACAGCTGTTGGCACGTGTGGACCA
TCTAATGAAAGTGGGAAAGAATAACTTCAGACCACGCCCAAGGTGGAATCCAGTGTTGTAAGGATAGAA
CCTAAGAATCCACCACCACCATCAATTTTCAGGAATGGGATGGTCTAGTAAGGATAACCTTTGTTAGGA
AAAACAAGACACTCTCTGCTGCATTTAAATCAAGTGCAGTGAACAACCTTTGAAAAAACTACAGAAT
TCACTGTTCACTCCATAATATTATAATACCAGAAGATTTCAGCATAGCAGATAAAATACAGCAAACTCTA
ACCAGCACAGGTTTTAGTGACAAACGGGCCCGTTCCATGGACATAGATGACTTCATCAGATTGCTACATG
GATTCAACGCAGAAGGTATTCATTTTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC202326 protein sequence Red=Cloning site Green=Tags(s) MPKVKSGAIGRRRGRQEQRRELKSAGGLMFNTGIGQHILKNPLIINSIIDKAALRPTDVLVLEVGPGTGNM TVKLLKAKKVACELDPRLVAELHKRVQGTPVASKLQVLVGDVLKTDLPFFDTCVANLPYQISSPFVFK LLLHRPFFRCAILMFQREFALRLVAKPGDKLYCRLSINTQLLARVDHLMKVGKNNFRPPPKVESSVRIE PKNPPPPINFQEWDLVRITFVRKNKTLAAFKSSAVQQLLEKNYRIHCSVHNIIPEDFSIADKIQQIL TSTGFSDKRARSMDIDDFIRLLHGFNAEGIHFS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_014473
ORF Size:	939 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_014473.2</u>
RefSeq Size:	1557 bp
RefSeq ORF:	942 bp
Locus ID:	27292
UniProt ID:	<u>Q9UNQ2</u>
Cytogenetics:	5q12.1
MW:	35.2 kDa