

Product datasheet for RC202328

PBLD (NM_022129) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PBLD (NM_022129) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PBLD
Synonyms:	HEL-S-306; MAWBP; MAWDBP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC202328 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**C

ATGAAGCTTCCTATTTTCATAGCAGATGCATTCACAGCAAGAGCATTTCTGTTGGAATCCTGCTGCTGTTT
 GCCTCCTAGAAAATGAATTGGATGAAGACATGCATCAGAAAATTGCAAGGGAGATGAACCTCTCTGAAAC
 TGCTTTTATCCGAAAACGCACCCGACAGACAACCTTGCACAAAGTTCCTGCTTTGGACTGAGATGGTTT
 ACACCAGCGAGTGAGTCCCACTCTGTGGCCATGCCACCCTGGCTTCTGCAGCTGTGCTGTTTCACAAAA
 TAAAAACATGAATAGCACGCTCACGTTTGTCACTCTGAGTGGAGAACTAAGGGCCAGACGAGCAGAGGA
 CGGCATCGTCCTGGACTTGCCTCTTTATCCAGCCACCCCGAGGACTTCCATGAAGTAGAGGACTTGATA
 AAGACTGCCATAGGCAACACACTGGTCCAGGACATCTGTTATTCTCCAGATACCCAAAAGCTCCTCGTCC
 GCCTCAGTGACGTTTACAACAGGTCGTTTCTGGAGAACCTGAAAGTGAACACGGAGAATCTGCTGCAAGT
 TGAAAACACAGGGAAGGTGAAAGGGCTTATTCTTACCCTTAAAGGAGAGCCTGGTGGGCAGACCAAGCA
 TTTGACTTTTACTCAAGATATTTTGCACCGTGGGTTGGTGTGGCTGAAGACCCAGTGACAGGGTCTGCAC
 ACGCTGTTCTCAGCAGCTACTGGTCCAGCATCTGGGGAAGAAAGAAATGCATGCTTTTCAGTGTCCCA
 CCGAGGAGGAGAGCTGGGAATTTCCCTTCGTCCAGACGGAAGGGTTGACATTAGAGGAGGTGCAGCTGTT
 GTTTTAGAGGGCACACTGACAGCC

ACGCGTACGCGGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >RC202328 protein sequence
 Red=Cloning site Green=Tags(s)


MKLPIFIADAFTARAFRGNPAAVCLLENELDEDMHQKIAREMNLSETAFIRKLHPTDNFAQSSCFGLRWF
 TPASEVPLCGHATLASAAVLHFKIKNMNSTLTFVTLSGELRARRAEDGIVLDLPLYPHPQDFHEVEDLI
 KTAIGNTLVQDICYSPDTQKLLVRLSDVYNRSFLENLKVNTENLLQVENTGKVKGLILTLKGEPGGQTQA
 FDFYSRYFAPWVGVAEDPVTGSAHAVLSSYWSQHLGKKEMHAFQCSHRGGELGISLRPDGRVDIRGGAAV
 VLEGLTA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6308_e04.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_022129

ORF Size: 864 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_022129.2](#)

RefSeq Size: 2587 bp

RefSeq ORF: 867 bp

Locus ID: 64081


UniProt ID: [P30039](#)

Cytogenetics: 10q21.3


Domains: PhzC-PhzF

MW: 31.8 kDa


Product images:




Circular map for RC202328



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY PBLD (Cat# RC202328, Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PBLD (Cat# [TA501994]). Positive lysates [LY411757] (100ug) and [LC411757] (20ug) can be purchased separately from OriGene.



Western blot validation of overexpression lysate (Cat# [LY411757]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC202328 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified PBLD protein (Cat# [TP302328]). The protein was produced from HEK293T cells transfected with PBLD cDNA clone (Cat# RC202328) using MegaTran 2.0 (Cat# [TT210002]).