

## Product datasheet for **RC202341L4V**

### ARMC8 (NM\_014154) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ARMC8 (NM_014154) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ARMC8
Synonyms:	GID5; HSPC056; S863-2; VID28
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_014154
ORF Size:	1155 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202341).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_014154.2</a>
RefSeq Size:	3236 bp
RefSeq ORF:	1158 bp
Locus ID:	25852
UniProt ID:	<a href="#">Q8IUR7</a>
Cytogenetics:	3q22.3
Domains:	Armadillo_seg
MW:	43 kDa



[View online »](#)

**Gene Summary:**

Component of the CTLH E3 ubiquitin-protein ligase complex that selectively accepts ubiquitin from UBE2H and mediates ubiquitination and subsequent proteasomal degradation of the transcription factor HBP1.[UniProtKB/Swiss-Prot Function]