

Product datasheet for **RC202346L3V**

LYRM4 (NM_020408) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	LYRM4 (NM_020408) Human Tagged ORF Clone Lentiviral Particle
Symbol:	LYRM4
Synonyms:	C6orf149; CGI-203; COXPD19; ISD11
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_020408
ORF Size:	273 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202346).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_020408.3
RefSeq Size:	1514 bp
RefSeq ORF:	276 bp
Locus ID:	57128
UniProt ID:	Q9HD34
Cytogenetics:	6p25.1
MW:	10.7 kDa


[View online »](#)

Gene Summary:

The protein encoded by this gene is found in both mitochondria and the nucleus, where it binds cysteine desulfurase and helps free inorganic sulfur for Fe/S clusters. Disruption of this gene negatively impacts mitochondrial and cytosolic iron homeostasis. [provided by RefSeq, Sep 2016]