

Product datasheet for **RC202359A1V**

Human TRA2A (NM_013293) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human TRA2A (NM_013293) AAV Particle
Tag:	Myc-DDK
Symbol:	TRA2A
Synonyms:	AWMS1; HSU53209
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC202359 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAGTGATGTGGAGGAAAACAACCTTCGAGGGCAGAGAGTCTCGCTCTCAGTCAAAATCTCCAACGGGAA
CTCCTGCTCGTGTAATTCGGAGAGCAGGTCAAGATCTCGTAGTCCATCAAGGGTTTCCAACACTCTGA
ATCCATTCTCGATCAAGATCAAAATCCAGGTCAAGGAGACATTCTCATAGACGTTACACTCGA
TCCAGATCCCACTCTCACTCTCATAGGAGACGATCTCGAAGTAGATCATATACACCAGAATACCGCGCGC
GAAGGAGCCGAAGCCATTCTCCAATGTCTAACCAGGAGAAGACATACTGGCAGCAGGGCAAATCCAGATCC
CAACACTTGCCTTGGAGTGTTTGGCCTCAGTTTGTACACAACAGAGAGGGATCTTCGTGAAGTATTTTCT
CGATATGGACCATTGAGTGGTGTCAATGTGTTTATGATCAGCGAACTGGGCGATCTCGAGGATTTGCTT
TTGTGATTTTGGAGAGATAGATGACTCAAAGGAGGCTATGGAAGGGCAAATGGAATGGAGCTGGATGG
TAGAAGAATTCGGGTGGATTATTCTATAACCAAGAGAGCGCACACCAACACCAGGCATCTACATGGGC
AGACCAACTCATAGTGGTGGGGTGGTGGAGGAGGCGCGCGGTGGAGGTGGAGGTGGTGGCAGACGTC
GAGATTCTTACTATGATAGAGGATATGATCGTGGGTATGACAGATATGAAGACTATGATTACCGATACAG
AAGACGATCACCTTCTCCTTATTATAGTCGATATAGATCAGATCAAGATCTCGTTCCTACAGCCCAAGA
CGCTAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC202359 protein sequence Red=Cloning site Green=Tags(s) MSDVEENNFEGRSRSQSKSPTGTPARVKSESRSRSGSRSPSRVSKHSESHSRSRSKSRSRRRHSHRRYTR SRSHSHSHRRRSRSRSTPEYRRRRSRSHSPMSNRRRHTGSRANPDNTCLGVFGLSLYTTERDLREVFS RYGPLSGVNVVYDQRTGRSRGFAFVYFERIDDSKEAMERANGMELDGRRIRVDYSITKRAHTPTPGIYMG RPTHSGGGGGGGGGGGGGGGRRRDSYYDRGYDRGYDRYEDYDYYRRRSPSPYYSRYRSRSRSRYSYSPR RY TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_013293
ORF Size:	846 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_013293.3</u>
RefSeq Size:	1866 bp
RefSeq ORF:	849 bp
Locus ID:	29896
UniProt ID:	<u>Q13595</u>
Cytogenetics:	7p15.3
MW:	32.7 kDa