

Product datasheet for **RC202359L3V**

TRA2A (NM_013293) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TRA2A (NM_013293) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TRA2A
Synonyms:	AWMS1; HSU53209
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_013293
ORF Size:	846 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202359).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013293.3
RefSeq Size:	1866 bp
RefSeq ORF:	849 bp
Locus ID:	29896
UniProt ID:	Q13595
Cytogenetics:	7p15.3
Domains:	RRM
Protein Pathways:	Spliceosome



[View online »](#)

MW: 32.7 kDa

Gene Summary: This gene is a member of the transformer 2 homolog family and encodes a protein with several RRM (RNA recognition motif) domains. This phosphorylated nuclear protein binds to specific RNA sequences and plays a role in the regulation of pre-mRNA splicing. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Sep 2013]