

Product datasheet for **RC202361L3V**

SLBP (NM_006527) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | SLBP (NM_006527) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | SLBP |
| Synonyms: | HBP |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_006527 |
| ORF Size: | 810 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC202361). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_006527.2 |
| RefSeq Size: | 1743 bp |
| RefSeq ORF: | 813 bp |
| Locus ID: | 7884 |
| UniProt ID: | Q14493 |
| Cytogenetics: | 4p16.3 |
| MW: | 31.3 kDa |


[View online »](#)

Gene Summary:

This gene encodes a protein that binds to the stem-loop structure in replication-dependent histone mRNAs. Histone mRNAs do not contain introns or polyadenylation signals, and are processed by endonucleolytic cleavage. The stem-loop structure is essential for efficient processing but this structure also controls the transport, translation and stability of histone mRNAs. Expression of the protein is regulated during the cell cycle, increasing more than 10-fold during the latter part of G1. [provided by RefSeq, Jul 2008]