

Product datasheet for **RC202364A1V**

Human DPF2 (NM_006268) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human DPF2 (NM_006268) AAV Particle
Tag:	Myc-DDK
Symbol:	DPF2
Synonyms:	CSS7; REQ; ubi-d4; UBID4
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC202364 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCGGCTGTGGTGGAGAATGTAGTGAAGCTCCTTGGGAGCAGTACTACAAAGATGCCATGGAGCAGT
GCCACAATTACAATGCTCGCCTCTGTGCTGAGCGCAGCGTGCGCCTGCCTTTCTTGACTCACAGACCGG
AGTAGCCAGAGCAATTGTTACATCTGGATGGAAAAGCGACACCGGGTCCAGGATTGGCCTCCGGACAG
CTGTACTCCTACCCTGCCCGGCGCTGGCGGAAAAAGCGCGAGCCATCCCCTGAGGATCCACGACTTT
CCTTCCCATCTATTAAGCCAGACACAGACCAGACCTGAAGAAGGAGGGGCTGATCTCTCAGGATGGCAG
TAGTTTAGAGGCTCTGTTGCGCACTGACCCCTGGAGAAGCGAGGTGCCCGGATCCCCGAGTTGATGAT
GACAGCCTGGGCGAGTTTCTGTGACCAACAGTCGAGCGCGAAAAGCGGATCCTAGAACCAGATGACTTCC
TGGATGACCTCGATGATGAAGACTATGAAGAAGATACTCCAAGCGTCGGGAAAGGGAAATCCAAGGG
TAAGGGTGTGGGAGTGCCCGTAAGAAGCTGGATGCTTCCATCCTGGAGGACCGGATAAGCCCTATGCC
TGTGACATTTGTGAAAACGTTACAAGAACCGACAGGCTCAGTTACCACTATGCCACTCCCCTTGG
CTGAGGAGGAGGCGAGGACAAGGAAGACTCTCAACCACTCCTGTTTCCAGAGGTCTGAGGAGCA
GAAATCAAAAAGGGTCTGATGGATTGGCCTTGCCCACTACTGTGACTTCTGCCTGGGGGACTCA
AAGATTAACAAGAAGACGGGACAACCCGAGGAGCTGGTGTCTGCTGACTGTGCCGCTCAGGCGATC
CATCTTGCCTCCAATTTACCCCGTGATGATGGCGGAGTGAAGACATACCGCTGGCAGTGCATCGAGTG
CAAATGTTGCAATATCTGCGGCACCTCCGAGAATGACGACCAAGTTGCTTCTGTGATGACTGCGATCGT
GGCTACCACATGTAAGTCTCACCCCGTCCATGTCTGAGCCCCCTGAAGGAAGTTGGAGCTGCCACCTGT
GTCTGGACCTGTTGAAAGAGAAAGCTTCCATCTACCAGAACCAAGTCTCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC202364 protein sequence Red=Cloning site Green=Tags(s) MAAVVENVVKLLGEQYYKDAMEQCHNYNARLCAERSVRLPFLDSQTGVAQSNCYIWMEKRHRGPGLASGQ LYSYPARRWKRRRAHPPEDPRLSFPSIKPDTDTLKKLEGLISQDGSSEALLRTDPLEKRGAPDPRVDD DSLGEFPVTNSRARKRILEPDDFLDDLDEDEYEDTPKRRGKSKGKGVGSARKKLDASILEDKPYA CDICGKRYKNRPGLSYHYAHSHLAEEEGEDKEDSQPPTPVSRSEEQKSKKGPDLALPNNYCDFCLGDS KINKKTGQPEELVSCSDCGRSGHPSCSQFTPVMMAAVKTYRWQCIECKCCNICGTSENDQQLFCDDCDR GYHMYCLTPSMSEPPEGSWSCHLCLDLLKEKASIYQNQNSS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_006268
ORF Size:	1173 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_006268.3
RefSeq Size:	2545 bp
RefSeq ORF:	1176 bp
Locus ID:	5977
UniProt ID:	Q92785
Cytogenetics:	11q13.1
MW:	44.2 kDa