

Product datasheet for **RC202386L3V**

PHF11 (NM_001040443) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | PHF11 (NM_001040443) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | PHF11 |
| Synonyms: | APY; BCAP; IGEL; IGER; IGHER; NY-REN-34; NYREN34 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_001040443 |
| ORF Size: | 876 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC202386). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001040443.1 |
| RefSeq Size: | 1649 bp |
| RefSeq ORF: | 996 bp |
| Locus ID: | 51131 |
| UniProt ID: | Q9UIL8 |
| Cytogenetics: | 13q14.2 |
| Protein Families: | Druggable Genome, Transcription Factors |
| MW: | 33.5 kDa |



[View online »](#)

Gene Summary:

This gene encodes a protein containing a PHD (plant homeodomain) type zinc finger. This gene has been identified in some studies as a candidate gene for asthma. Naturally-occurring readthrough transcription may occur from the upstream SETDB2 (SET domain bifurcated 2) gene to this locus. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2016]