

Product datasheet for **RC202389L4V**

Frequenin (NCS1) (NM_014286) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Frequenin (NCS1) (NM_014286) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Frequenin
Synonyms:	FLUP; FREQ
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_014286
ORF Size:	570 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202389).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_014286.2 , NP_055101.2
RefSeq Size:	5009 bp
RefSeq ORF:	573 bp
Locus ID:	23413
UniProt ID:	P62166
Cytogenetics:	9q34.11
Domains:	EFh
Protein Families:	Druggable Genome


[View online »](#)

MW: 21.9 kDa

Gene Summary: This gene is a member of the neuronal calcium sensor gene family, which encode calcium-binding proteins expressed predominantly in neurons. The protein encoded by this gene regulates G protein-coupled receptor phosphorylation in a calcium-dependent manner and can substitute for calmodulin. The protein is associated with secretory granules and modulates synaptic transmission and synaptic plasticity. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]