

Product datasheet for **RC202394L2V**

SYNJ2BP (NM_018373) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SYNJ2BP (NM_018373) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SYNJ2BP
Synonyms:	ARIP2; OMP25
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_018373
ORF Size:	435 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202394).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018373.1
RefSeq Size:	7074 bp
RefSeq ORF:	438 bp
Locus ID:	55333
UniProt ID:	P57105
Cytogenetics:	14q24.2
Protein Families:	Transmembrane
MW:	15.9 kDa



[View online »](#)

Gene Summary:

Regulates endocytosis of activin type 2 receptor kinases through the Ral/RALBP1-dependent pathway and may be involved in suppression of activin-induced signal transduction.
[UniProtKB/Swiss-Prot Function]