

Product datasheet for **RC202398A1V**

Human GOLPH3L (NM_018178) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human GOLPH3L (NM_018178) AAV Particle
Tag:	Myc-DDK
Symbol:	GOLPH3L
Synonyms:	GPP34R
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC202398 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGACCACTTTAACTACCGGGCCCGTCGCACTGAAATAAGCAAGAACTCTGAAAAGAAGATGGAAAGTG
AGGAAGACAGTAATTGGGAGAAAAGTCCAGACAATGAAGATTCTGGAGACTCTAAGGATATCCGCCTTAC
TCTTATGGAAGAAGTATTGCTTCTGGGACTAAAAGATAAAGAGGGGTACACATCTTTCTGGAATGACTGC
ATATCATCAGGCCTGCGAGGGGGCATCCTGATAGAGCTGGCCATGCGGGGTGCAATCTATCTGGAACCCC
CGACCATGCGTAAGAAGCGACTACTAGACAGAAAGGTACTGCTAAAGTCAGACAGCCCAACAGGTGATGT
TTTACTGGATGAACTCTGAAACACATCAAAGCAACTGAACCCACAGAACTGTCCAAACATGGATAGAG
CTACTCACTGGTGAGACCTGGAACCCCTTCAAATTACAGTACCAGCTGAGAAATGTACGAGAGCGCATCG
CAAAGAACCTAGTAGAGAAAGGTATTCTAACCCTGAGAAGCAGAATTTCTGCTATTTGACATGACTAC
TCATCCAGTGACCAATACAACAGAGAAACAGCGACTAGTGAAAAAACTTCAAGATAGTGTACTAGAGCGG
TGGGTAAATGACCCTCAGCGTATGGACAAGCGAACTAGCACTCCTGGTGCTAGCCCACTCCTCTGATG
TGCTAGAGAATGTCTTCTCTCTGACAGATGACAAGTATGATGTGGCAATGAATCGAGCCAAGGACTT
AGTAGAACTGGACCCTGAAGTGAAGGGACAAAGCCTAGTGCCACAGAAATGATCTGGGCTGTGCTGGCA
GCCTTCAATAAATCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC202398 protein sequence Red=Cloning site Green=Tags(s) MTTLTHRARRTEISKNSEKKMESEEDSNWEKSPDNEDSGDSKDIRLTLMEEVLLLGLKDKEGYTSFWNDC ISSGLRGGILIELAMRGRIYLEPPTMRKKRLDRKVLLKSDSPTGDVLLDETLKHIKATEPTETVQTWIE LLTGETWNPFLQYQLRNVREIRIAKNLVEKGILTTEKQNFLLFDMTHPVTNTTEKQRLVKKLQDSVLER WYNDPQRMDKRTLALLVLAHSSDVLENVSSLTDDKYDVAMNRAKDLVELDPEVEGTPKSATEMIWAVLA AFNKS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_018178
ORF Size:	855 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_018178.3</u>
RefSeq Size:	3208 bp
RefSeq ORF:	858 bp
Locus ID:	55204
UniProt ID:	<u>Q9H4A5</u>
Cytogenetics:	1q21.3
MW:	32.8 kDa