

## Product datasheet for **RC202410L3V**

### APMAP (NM\_020531) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	APMAP (NM_020531) Human Tagged ORF Clone Lentiviral Particle
Symbol:	APMAP
Synonyms:	BSCv; C20orf3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_020531
ORF Size:	1248 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202410).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_020531.2</a>
RefSeq Size:	2251 bp
RefSeq ORF:	1251 bp
Locus ID:	57136
UniProt ID:	<a href="#">Q9HDC9</a>
Cytogenetics:	20p11.21
Domains:	Str_synth
Protein Families:	Transmembrane



[View online »](#)

**MW:** 46.5 kDa

**Gene Summary:** Exhibits strong arylesterase activity with beta-naphthyl acetate and phenyl acetate. May play a role in adipocyte differentiation.[UniProtKB/Swiss-Prot Function]