

## Product datasheet for **RC202428L3V**

### SUV39H1 (NM\_003173) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:             | Lentiviral Particles                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Product Name:             | SUV39H1 (NM_003173) Human Tagged ORF Clone Lentiviral Particle                                                                                                                                                                                                                                                                                                                                                                                     |
| Symbol:                   | SUV39H1                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Synonyms:                 | H3-K9-HMTase 1; KMT1A; MG44; SUV39H                                                                                                                                                                                                                                                                                                                                                                                                                |
| Mammalian Cell Selection: | Puromycin                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Vector:                   | pLenti-C-Myc-DDK-P2A-Puro (PS100092)                                                                                                                                                                                                                                                                                                                                                                                                               |
| Tag:                      | Myc-DDK                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| ACCN:                     | NM_003173                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| ORF Size:                 | 1236 bp                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC202428).                                                                                                                                                                                                                                                                                                                                                                                     |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.                                                                                                                                                                                                                                                                                                               |
| RefSeq:                   | <a href="#">NM_003173.2</a>                                                                                                                                                                                                                                                                                                                                                                                                                        |
| RefSeq Size:              | 2752 bp                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| RefSeq ORF:               | 1239 bp                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Locus ID:                 | 6839                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| UniProt ID:               | <a href="#">O43463</a>                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Cytogenetics:             | Xp11.23                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Domains:                  | CHROMO, SET, PreSET, PostSET, Pre-SET                                                                                                                                                                                                                                                                                                                                                                                                              |
| Protein Families:         | Druggable Genome                                                                                                                                                                                                                                                                                                                                                                                                                                   |



[View online »](#)

**Protein Pathways:** Lysine degradation

**MW:** 47.9 kDa

**Gene Summary:** This gene encodes an evolutionarily-conserved protein containing an N-terminal chromodomain and a C-terminal SET domain. The encoded protein is a histone methyltransferase that trimethylates lysine 9 of histone H3, which results in transcriptional gene silencing. Loss of function of this gene disrupts heterochromatin formation and may cause chromosome instability. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2013]