

Product datasheet for **RC202429L4V**

SCGF (CLEC11A) (NM_002975) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SCGF (CLEC11A) (NM_002975) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SCGF
Synonyms:	CLECSF3; LSLCL; P47; SCGF
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_002975
ORF Size:	969 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202429).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002975.2
RefSeq Size:	1423 bp
RefSeq ORF:	972 bp
Locus ID:	6320
UniProt ID:	Q9Y240
Cytogenetics:	19q13.33
Protein Families:	Druggable Genome, Secreted Protein
MW:	35.7 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of the C-type lectin superfamily. The encoded protein is a secreted sulfated glycoprotein and functions as a growth factor for primitive hematopoietic progenitor cells. An alternative splice variant has been described but its biological nature has not been determined. [provided by RefSeq, Jul 2008]