

Product datasheet for **RC202441L2V**

CAMKV (NM_024046) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CAMKV (NM_024046) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CAMKV
Synonyms:	1G5; VACAMKL
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_024046
ORF Size:	1503 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202441).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024046.2
RefSeq Size:	3028 bp
RefSeq ORF:	1506 bp
Locus ID:	79012
UniProt ID:	Q8NCB2
Cytogenetics:	3p21.31
Domains:	ptkase, TyrKc, S_TKc
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)

MW: 54.4 kDa

Gene Summary: Does not appear to have detectable kinase activity.[UniProtKB/Swiss-Prot Function]