

Product datasheet for **RC202447A1V**

Human IL1RA (IL1RN) (NM_173843) AAV Particle

Product data:

Product Type: AAV Particles

Product Name: Human IL1RA (IL1RN) (NM_173843) AAV Particle

Tag: Myc-DDK

Symbol: IL1RA

Synonyms: DIRA; ICIL-1RA; IL-1ra; IL-1ra3; IL-1RN; IL1F3; IL1RA; IRAP; MVCD4

Mammalian Cell Selection: None

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC202447 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCTTTAGAGACGATCTGCCGACCCTCTGGGAGAAAATCCAGCAAGATGCAAGCCTTCAGAATCTGGG
ATGTTAACCAGAAGACCTTCTATCTGAGGAACAACCACTAGTTGCCGGATACTTGCAAGGACCAATGT
CAATTTAGAAGAAAAGATAGATGTGGTACCCATTGAGCCTCATGCTCTGTTCTTGGGAATCCATGGAGGG
AAGATGTGCCTGTCCTGTGTCAAGTCTGGTGATGAGACCAGACTCCAGCTGGAGGCAGTTAACATCACTG
ACCTGAGCGAGAACAGAAAGCAGGACAAGCGCTTCGCCTTCATCCGCTCAGACAGTGGCCCCACCACAG
TTTTGAGTCTGCCGCTGCCCGGTTGGTTCTCTGCACAGCGATGGAAGCTGACCAGCCCGTCAGCCTC
ACCAATATGCCTGACGAAGCGTCATGGTCACCAAATTCTACTTCCAGGAGGACGAG

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC202447 protein sequence
Red=Cloning site Green=Tags(s)

MALETICRPSGRKSSKMQAFLRIWDVNQKTFYLRNNQLVAGYLQGPVNLEEKIDVVPPIEPHALFLGIHGG
KMCLSCVKSGDETRLQLEAVNITDLSENKQDKRFAFIRSDSGPTTSFESAACPGWFLCTAMEADQPVSL
TNMPDEGVMVTKFYFQEDE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human

Serotype: AAV-2



[View online »](#)

ACCN:	NM_173843
ORF Size:	477 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_173843.1</u>
RefSeq Size:	1973 bp
RefSeq ORF:	432 bp
Locus ID:	3557
UniProt ID:	<u>P18510</u>
Cytogenetics:	2q14.1
MW:	17.9 kDa