

Product datasheet for **RC202449A1V**

Human beta Synuclein (SNCB) (NM_003085) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human beta Synuclein (SNCB) (NM_003085) AAV Particle
Tag: Myc-DDK
Symbol: beta Synuclein
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC202449 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGACGTGTTTCATGAAGGGCCTGTCCATGGCCAAGGAGGGCGTCTGGCAGCCGCGGAGAAAACCAAGC
 AGGGGGTCACCGAGGCGCGGAGAAGACCAAGGAGGGCGTCTCTACGTCGGAAGCAAGACCCGAGAAGG
 TGTGGTACAAGGTGTGGCTTCAGTGGCTGAAAAACCAAGGAACAGGCCTCACATCTGGGAGGAGCTGTG
 TTCTCTGGGGCAGGGAACATCGCAGCAGCCACAGGACTGGTGAAGAGGGAGGAATTCCTACTGATCTGA
 AGCCAGAGGAAGTGGCCAGGAAGCTGCTGAAGAACCCTGATTGAGCCCTGATGGAGCCAGAAGGGGA
 GAGTTATGAGGACCCACCCAGGAGGAATATCAGGAGTATGAGCCAGAGGCG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC202449 protein sequence
 Red=Cloning site Green=Tags(s)

MDVFMKGLSMAKEGVVAAAETKQGVTEAAETKKEGVLVYVSKTREGVVQGVASVAEKTKEQASHLGGAV
 FSGAGNIAAATGLVKREEFPDLDKPEEVAQEAEEPLIEPLMEPEGESYEDPPQEYQEYEPEA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_003085
ORF Size: 402 bp



[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_003085.3</u>
RefSeq Size:	1310 bp
RefSeq ORF:	405 bp
Locus ID:	6620
UniProt ID:	<u>Q16143</u>
Cytogenetics:	5q35.2
MW:	14.3 kDa