

Product datasheet for **RC202457A1V**

Human CC130 (CCDC130) (NM_030818) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human CC130 (CCDC130) (NM_030818) AAV Particle
Tag:	Myc-DDK
Symbol:	CC130
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC202457 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGGTGAAAGGAAAGGGGTCAACAAGTACTATCCTCCAGACTTCAACCCTGAGAAGCATGGCTCTCTCA
ACCGATACCACAACAGCCACCCGCTTCGGGAGCGGGCTCGGAAGCTGTACAGGGCATCCTCATCATCCG
ATTCGAAATGCCATATAACATCTGGTGCATGGCTGCAAGAACCACATCGGCATGGGTGTTCTGTTACAAT
GCAGAAAAGAAGAAGGTGGGCAATTACTACACAACCCCGATCTACAGGTTCCGGATGAAATGCCACCTCT
GTGTCAACTACATCGAGATGCAGACGGACCCCGCAACTGCGACTACGTGATCGTGAGTGGCGCCACGCG
CAAGGAGGAGCGCTGGGACATGGCGGACAATGAGCAGGTGCTGACCACAGAGCATGAGAAGAAGCAGAAG
CTGGAGACGGACGCCATGTTCCGGCTGGAGCATGGCGAGGCCGACCGCAGCACACTCAAGAAGGCGCTGC
CCACACTGAGCCACATCCAGGAGGCCAGAGCGCTGGAAGGACGACTTCGCCCTCAACAGCATGCTGCG
GAGAAGGTTCCGGGAAAAGAAAAAGCCATCCAGGAGGAGGAGAGAGACCAGGCCTTGCAAGGCCAAG
GCGAGCCTGACCATCCCGCTGGTGCCGAGACGGAAGTGACCGCAAGCTGGCGGCTCTGCTGAAGTTCC
ACACCTGGACTCCTACGAGGACAAGCAGAACTCAAGCGGACCGAGATCATCAGCCGCTCCTGGTTCCC
CTTGCCCCCGGATCCGCCTCCAGCAGCAAGGTCAGCGGCTCCTGAAGAAGCTGGCACAGAGCCGAGAG
ACCGCGCTTGCCACCTCCCCATCACCGTCGGGACCTGGGCATCGTGGCGGAGGTCTCGGACGTCC
CGGAGAGCCCCCAGCATGCGGCCGACACCCCCAAGTCTGGGGAACCGCGGTACCAGAGGAGGTGCCCA
GGACCGGCCATGTCCCCCGGAGACTGTCTCCGAAACAAGTGAACCCCAAGTGCAGCAGCCGAGG
GGGAGGAAGGAGCGGTGAGGACAAGCCCTGTGCGCAGCAGGCTCTCCAGGAGGCAGCTGACACCC
CCGACACGCGGACCCCTGCAGTCTCGGCTCCTCCCTCGTGGCGGACTACTCCGACTCGGAGAGTGAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC202457 protein sequence Red=Cloning site Green=Tags(s) MGERKGVNKYYPPDFNPEKHGSLNRYHNSHPLRERARKLSQGILIIIRFEMPYNIWCDGCKNHIGMGVRYN AEKKKVGNYTTPIYRFRMKCHLCVNYIEMQTDPANCDYVIVSGAQRKEERWDMADNEQVL TTEHEKKQK LETDA MFRLEHGEADRSTLKKALPTLSHIQEAQSAWKDDFALNSMLRRRFREKKKAIQEEEEERDQALQAK ASLTIP LVPETEDDRKLAALLKFHTLDSYEDKQKLK RTEIISRWFPSAPGSASSSKVSGVLKKLAQSRR TALATSPITVGD LGIVRRRSRDVPES PQHAADTPKSGEPRVP EEAQDRPMSPGDCPPETTETPKCSPR GQEGSRQDKPLSPAGSSQEAADTPDTRHPCSLGSSLVADYSDSESE TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_030818
ORF Size:	1188 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_030818.2</u>
RefSeq Size:	1679 bp
RefSeq ORF:	1191 bp
Locus ID:	81576
UniProt ID:	<u>P13994</u>
Cytogenetics:	19p13.13
MW:	44.8 kDa