

# Product datasheet for RC202460L2V

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## MMP1 (NM\_002421) Human Tagged ORF Clone Lentiviral Particle

#### **Product data:**

Product Type: Lentiviral Particles

**Product Name:** MMP1 (NM\_002421) Human Tagged ORF Clone Lentiviral Particle

Symbol: MMP1

Synonyms: CLG; CLGN

Mammalian Cell

Selection:

None

**Vector:** pLenti-C-mGFP (PS100071)

Tag: mGFP

**ACCN:** NM\_002421 **ORF Size:** 1407 bp

**ORF Nucleotide** 

The ORF insert of this clone is exactly the same as(RC202460).

OTI Disclaimer:

Sequence:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**RefSeg:** NM 002421.2, NP 002412.1

 RefSeq Size:
 2081 bp

 RefSeq ORF:
 1410 bp

 Locus ID:
 4312

 UniProt ID:
 P03956

 Cytogenetics:
 11q22.2

**Domains:** hemopexin, Peptidase\_M10, ZnMc

**Protein Families:** Druggable Genome, Protease, Secreted Protein





### MMP1 (NM\_002421) Human Tagged ORF Clone Lentiviral Particle - RC202460L2V

**Protein Pathways:** Bladder cancer, Pathways in cancer, PPAR signaling pathway

**MW:** 54 kDa

**Gene Summary:** This gene encodes a member of the peptidase M10 family of matrix metalloproteinases

(MMPs). Proteins in this family are involved in the breakdown of extracellular matrix in normal physiological processes, such as embryonic development, reproduction, and tissue remodeling, as well as in disease processes, such as arthritis and metastasis. The encoded preproprotein is proteolytically processed to generate the mature protease. This secreted protease breaks down the interstitial collagens, including types I, II, and III. The gene is part of a cluster of MMP genes on chromosome 11. Mutations in this gene are associated with chronic obstructive pulmonary disease (COPD). Alternative splicing results in multiple transcript variants, at least one of which encodes an isoform that is proteolytically processed.

[provided by RefSeq, Jan 2016]