

Product datasheet for RC202470L1

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Optineurin (OPTN) (NM_001008211) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: Optineurin (OPTN) (NM 001008211) Human Tagged Lenti ORF Clone

Tag: Myc-DDK
Symbol: Optineurin

Synonyms: ALS12; FIP2; GLC1E; HIP7; HYPL; NRP; TFIIIA-INTP

Mammalian Cell None

Selection:

Vector:pLenti-C-Myc-DDK (PS100064)E. coli Selection:Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC202470).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF

ACCN: NM_001008211

ORF Size: 1731 bp





Optineurin (OPTN) (NM_001008211) Human Tagged Lenti ORF Clone - RC202470L1

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 001008211.1</u>, <u>NP 001008212.1</u>

 RefSeq Size:
 3597 bp

 RefSeq ORF:
 1734 bp

 Locus ID:
 10133

 UniProt ID:
 Q96CV9

 Cytogenetics:
 10p13

Protein Families: Druggable Genome

MW: 65.9 kDa


Gene Summary: This gene encodes the coiled-coil containing protein optineurin. Optineurin may play a role in

normal-tension glaucoma and adult-onset primary open angle glaucoma. Optineurin interacts with adenovirus E3-14.7K protein and may utilize tumor necrosis factor-alpha or Fas-ligand pathways to mediate apoptosis, inflammation or vasoconstriction. Optineurin may also function in cellular morphogenesis and membrane trafficking, vesicle trafficking, and transcription activation through its interactions with the RAB8, huntingtin, and transcription factor IIIA proteins. Alternative splicing results in multiple transcript variants encoding the


same protein. [provided by RefSeq, Jul 2008]



Product images:



Circular map for RC202470L1



Double digestion of RC202470L1 using Sgfl and Mlul $\,$