

Product datasheet for **RC202473L3V**

UBE2D4 (NM_015983) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	UBE2D4 (NM_015983) Human Tagged ORF Clone Lentiviral Particle
Symbol:	UBE2D4
Synonyms:	HBUCE1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_015983
ORF Size:	441 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202473).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_015983.2
RefSeq Size:	1398 bp
RefSeq ORF:	444 bp
Locus ID:	51619
UniProt ID:	Q9Y2X8
Cytogenetics:	7p13
Domains:	UBCc
Protein Pathways:	Ubiquitin mediated proteolysis



[View online »](#)

MW: 16.5 kDa

Gene Summary: Accepts ubiquitin from the E1 complex and catalyzes its covalent attachment to other proteins. In vitro able to promote polyubiquitination using all 7 ubiquitin Lys residues, but may prefer 'Lys-11' and 'Lys-48'-linked polyubiquitination.[UniProtKB/Swiss-Prot Function]