

Product datasheet for **RC202486L2V**

EHD4 (NM_139265) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	EHD4 (NM_139265) Human Tagged ORF Clone Lentiviral Particle
Symbol:	EHD4
Synonyms:	PAST4
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_139265
ORF Size:	1623 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202486).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_139265.2 , NP_644670.1
RefSeq Size:	2904 bp
RefSeq ORF:	1626 bp
Locus ID:	30844
UniProt ID:	Q9H223
Cytogenetics:	15q15.1
Domains:	EH
Protein Pathways:	Endocytosis



[View online »](#)

MW: 61.2 kDa

Gene Summary: ATP- and membrane-binding protein that probably controls membrane reorganization/tubulation upon ATP hydrolysis. Plays a role in early endosomal transport. [UniProtKB/Swiss-Prot Function]