

Product datasheet for RC202487A1V

Human VPS25 (NM_032353) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human VPS25 (NM_032353) AAV Particle
Tag: Myc-DDK
Symbol: VPS25
Synonyms: DERP9; EAP20; FAP20
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC202487 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGGCGATGAGTTTCGAGTGGCCGTGGCAGTATCGCTTCCACCTTCTTTACGTTACAACCGAATGTGG
 ACATCTGGCAGAAGCAGCTGGCCGCTGGTGCTCGCTGGTCTGTCTTCTGCCGCTGCACAAACAGTC
 CAGCATGACGGTGATGGAAGCTCAGGAGAGCCGCTCTTCAACAACGTCAAGCTACAGCGAAAGCTTCT
 GTGGAGTCGATCCAGATTGTATTAGAGAACTGAGGAAGAAAGGAACCTCGAGTGGTTGGATAAGAGCA
 AGTCCAGCTTCTGATCATGTGGCGGAGGCCAGAAGAATGGGGAACTCATCTATCAGTGGGTTCCAG
 GAGTGGCCAGAACAACCTCCGCTTTACCCTGTATGAAGTGAATGAGGGAAGACACAGAGGATGAGGAG
 TTCCACGGGCTGGATGAAGCCACTCTACTGCGGCTCTGCAGGCCCTACAGCAGGAGCACAAGGCCGAGA
 TCATCACTGTACGCGATGGCCGAGGCGTCAAGTTCTTC

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC202487 protein sequence
 Red=Cloning site Green=Tags(s)

MAMSEFWPWQYRFPFFTLQPNVDTRQKQLAAWCSLVLSFCRLHKQSSMTVMEAQESPLFNNVKLQRKLP
 VESIQIVLEELRKKGNLEWLDKSKSSFLIMWRRPEEWGKLIYQWVSRSGQNNVFTLYELTNGEDTEDEE
 FHGLDEATLLRALQALQQEHKAEIITVSDGRGVKFF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2


[View online »](#)

ACCN:	NM_032353
ORF Size:	528 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_032353.2</u>
RefSeq Size:	1120 bp
RefSeq ORF:	531 bp
Locus ID:	84313
UniProt ID:	<u>Q9BRG1</u>
Cytogenetics:	17q21.2
MW:	20.7 kDa