

Product datasheet for **RC202496L1V**

ZNF43 (NM_003423) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF43 (NM_003423) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF43
Synonyms:	HTF6; KOX27; ZNF39L1
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_003423
ORF Size:	2427 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202496).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_003423.3
RefSeq Size:	5258 bp
RefSeq ORF:	2430 bp
Locus ID:	7594
UniProt ID:	P17038
Cytogenetics:	19p12
Domains:	KRAB, zf-C2H2
Protein Families:	Transcription Factors



[View online »](#)

MW: 94.1 kDa

Gene Summary: This gene belongs to the C2H2-type zinc finger gene family. The zinc finger proteins are involved in gene regulation and development, and are quite conserved throughout evolution. Like this gene product, a third of the zinc finger proteins containing C2H2 fingers also contain the KRAB domain, which has been found to be involved in protein-protein interactions. [provided by RefSeq, Jul 2008]