

Product datasheet for RC202497A1V

Human Uteroglobin (SCGB1A1) (NM_003357) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human Uteroglobin (SCGB1A1) (NM_003357) AAV Particle
Tag: Myc-DDK
Symbol: Uteroglobin
Synonyms: CC10; CC16; CCPBP; CCSP; UGB; UP-1; UP1
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC202497 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGAACTCGCTGTCAACCTCACCTGGTCACACTGGCTCTCTGCTGCAGCTCCGCTTCTGCAGAGATCT
 GCCCGAGCTTTTCAGCGTGTATCGAAACCCTCCTCATGGACACACCCTCCAGTTATGAGGCTGCCATGGA
 ACTTTTCAGCCCTGATCAAGACATGAGGGAGGCAGGGGCTCAGCTGAAGAAGCTGGTGGACACCTCCCC
 CAAAAGCCAGAGAAAGCATCATTAAAGCTCATGGAAAAATAGCCCAAAGCTCACTGTGTAAT

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC202497 protein sequence
 Red=Cloning site Green=Tags(s)

MKLAVTLTLVTALCCSSASAEICPSFQRVIETLLMDTPSSYEAMELFSPDQDMREAGQLKKLVDTLP
 QKPREIIKLMEKIAQSSLCN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_003357
ORF Size: 273 bp
Buffer: PBS with 0.001% Pluronic F68


[View online »](#)

Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_003357.3</u>
RefSeq Size:	452 bp
RefSeq ORF:	276 bp
Locus ID:	7356
UniProt ID:	<u>P11684</u>
Cytogenetics:	11q12.3
MW:	10 kDa