

## Product datasheet for **RC202497L1V**

### Uteroglobin (SCGB1A1) (NM\_003357) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Uteroglobin (SCGB1A1) (NM_003357) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Uteroglobin
Synonyms:	CC10; CC16; CCPBP; CCSP; UGB; UP-1; UP1
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_003357
ORF Size:	273 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202497).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_003357.3</a>
RefSeq Size:	452 bp
RefSeq ORF:	276 bp
Locus ID:	7356
UniProt ID:	<a href="#">P11684</a>
Cytogenetics:	11q12.3
Protein Families:	Druggable Genome, Secreted Protein
MW:	10 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of the secretoglobin family of small secreted proteins. The encoded protein has been implicated in numerous functions including anti-inflammation, inhibition of phospholipase A2 and the sequestering of hydrophobic ligands. Defects in this gene are associated with a susceptibility to asthma. [provided by RefSeq, May 2010]