

## Product datasheet for **RC202502L1V**

### Dopamine Receptor D5 (DRD5) (NM\_000798) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Dopamine Receptor D5 (DRD5) (NM_000798) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Dopamine Receptor D5
Synonyms:	DBDR; DRD1B; DRD1L2
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_000798
ORF Size:	1431 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202502).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_000798.3</a> , <a href="#">NP_000789.1</a>
RefSeq Size:	2398 bp
RefSeq ORF:	1434 bp
Locus ID:	1816
UniProt ID:	<a href="#">P21918</a>
Cytogenetics:	4p16.1
Protein Families:	Druggable Genome, GPCR, Transmembrane
Protein Pathways:	Calcium signaling pathway, Neuroactive ligand-receptor interaction


[View online »](#)

**MW:** 53 kDa

**Gene Summary:** This gene encodes the D5 subtype of the dopamine receptor. The D5 subtype is a G-protein coupled receptor which stimulates adenylyl cyclase. This receptor is expressed in neurons in the limbic regions of the brain. It has a 10-fold higher affinity for dopamine than the D1 subtype. Pseudogenes related to this gene reside on chromosomes 1 and 2. [provided by RefSeq, Jul 2008]