

## Product datasheet for RC202502L3

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Dopamine Receptor D5 (DRD5) (NM\_000798) Human Tagged Lenti ORF Clone

**Product data:** 

**Product Type:** Expression Plasmids

Product Name: Dopamine Receptor D5 (DRD5) (NM\_000798) Human Tagged Lenti ORF Clone

Tag: Myc-DDK

Symbol: Dopamine Receptor D5
Synonyms: DBDR; DRD1B; DRD1L2

Mammalian Cell Puromycin

Selection:

**Vector:** pLenti-C-Myc-DDK-P2A-Puro (PS100092)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC202502).

Sequence:

Sites: Sgfl-Mlul

Restriction Sites: Cloning Scheme:





<sup>\*</sup> The last codon before the Stop codon of the ORF.

**ACCN:** NM\_000798

ORF Size: 1431 bp



## Dopamine Receptor D5 (DRD5) (NM\_000798) Human Tagged Lenti ORF Clone - RC202502L3

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 000798.3</u>, <u>NP 000789.1</u>

 RefSeq Size:
 2398 bp

 RefSeq ORF:
 1434 bp

 Locus ID:
 1816

 UniProt ID:
 P21918

**Protein Families:** Druggable Genome, GPCR, Transmembrane

4p16.1

**Protein Pathways:** Calcium signaling pathway, Neuroactive ligand-receptor interaction

**MW:** 53 kDa

Cytogenetics:

**Gene Summary:** This gene encodes the D5 subtype of the dopamine receptor. The D5 subtype is a G-protein

coupled receptor which stimulates adenylyl cyclase. This receptor is expressed in neurons in the limbic regions of the brain. It has a 10-fold higher affinity for dopamine than the D1 subtype. Pseudogenes related to this gene reside on chromosomes 1 and 2. [provided by

RefSeq, Jul 2008]