

Product datasheet for **RC202516A1V**

Human UBE2H (NM_003344) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human UBE2H (NM_003344) AAV Particle
Tag: Myc-DDK
Symbol: UBE2H
Synonyms: E2-20K; GID3; UBC8; UBCH; UBCH2
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC202516 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTCATCTCCAGTCCGGGCAAGAGGCGGATGGACACGGACGTGGTCAAGCTCATCGAGAGTAAACATG
AGGTTACGATCCTGGGAGGACTTAATGAATTTGTAGTGAAGTTTTATGGACCACAAGGAACCATATGA
AGGCGGAGTATGGAAAGTTAGAGTGGACCTACCTGATAAATACCTTTCAAATCTCCATCTATAGGATTC
ATGAATAAAATTTCCATCCCAACATTGATGAAGCGTCAGGAAGTGTGTCTAGATGTAATTAATCAAA
CTTGGACAGCTCTCTATGATCTTACCAATATATTGAGTCCTTCCTGCCTCAGTTATTGGCTATCCTAA
CCCCATAGATCCTCTCAATGGTGACGCTGCAGCCATGTACCTCCACCGACCAGAAGAATACAAGCAGAAA
ATTAAAGAGTACATCCAGAAATACGCCACGGAGGAGCGCTGAAAGAACAGGAAGAGGGTACCGGGGACA
GCTCATCGGAGAGCTCTATGTCTGACTTTTCCGAAGATGAGGCCAGGATATGGAGTTG

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC202516 protein sequence
Red=Cloning site Green=Tags(s)

MSSSPGKRRMDTDVVKLIESKHEVTILGGLNEFVVKFYGPQGTPEGGVWKVRVDLPDKYFPKSPSIGF
MNKIFHPNIDEASGTVCLDVINQWTALYDLTNIFESFLPQLLAYPNPIDPLNGDAAAMYLHRPEEYKQK
IKEYIQKYATEEALKEQEETGDSSSESSMSDFSEDEAQMEL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2



[View online »](#)

ACCN:	NM_003344
ORF Size:	549 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_003344.2</u>
RefSeq Size:	5174 bp
RefSeq ORF:	552 bp
Locus ID:	7328
UniProt ID:	<u>P62256</u>
Cytogenetics:	7q32.2
MW:	20.7 kDa