

Product datasheet for **RC202517L4V**

HOOK1 (NM_015888) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	HOOK1 (NM_015888) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HOOK1
Synonyms:	HK1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_015888
ORF Size:	2184 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202517).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_015888.4
RefSeq Size:	5865 bp
RefSeq ORF:	2187 bp
Locus ID:	51361
UniProt ID:	Q9UJC3
Cytogenetics:	1p32.1
MW:	84.6 kDa

[View online »](#)

Gene Summary:

This gene encodes a member of the hook family of coiled-coil proteins, which bind to microtubules and organelles through their N- and C-terminal domains, respectively. The encoded protein localizes to discrete punctuate subcellular structures, and interacts with several members of the Rab GTPase family involved in endocytosis. It is thought to link endocytic membrane trafficking to the microtubule cytoskeleton. Several alternatively spliced transcript variants have been identified, but the full-length nature of some of these variants has not been determined. [provided by RefSeq, Jul 2008]