

## Product datasheet for RC202529L1

### NDUFAF7 (NM\_144736) Human Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	NDUFAF7 (NM_144736) Human Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	NDUFAF7
Synonyms:	C2orf56; MidA; PRO1853
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202529).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.


ACCN:	NM_144736
ORF Size:	1323 bp



[View online »](#)

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<a href="#">NM_144736.3</a>
<b>RefSeq Size:</b>	2221 bp
<b>RefSeq ORF:</b>	1326 bp
<b>Locus ID:</b>	55471
<b>UniProt ID:</b>	<a href="#">Q7L592</a>
<b>Cytogenetics:</b>	2p22.2
<b>Domains:</b>	DUF185
<b>MW:</b>	49.2 kDa
<b>Gene Summary:</b>	This gene encodes an assembly factor protein which helps in the assembly and stabilization of Complex I, a large multi-subunit enzyme in the mitochondrial respiratory chain. Complex I is involved in several physiological activities in the cell, including metabolite transport and ATP synthesis. The encoded protein is a methyltransferase which methylates Arg85 of a subunit of Complex I in the early stages of its assembly. A pseudogene related to this gene is located on chromosome 8. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2016]

## Product images:



Circular map for RC202529L1