

Product datasheet for **RC202529L1V**

NDUFAF7 (NM_144736) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	NDUFAF7 (NM_144736) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NDUFAF7
Synonyms:	C2orf56; MidA; PRO1853
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_144736
ORF Size:	1323 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202529).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_144736.3
RefSeq Size:	2221 bp
RefSeq ORF:	1326 bp
Locus ID:	55471
UniProt ID:	Q7L592
Cytogenetics:	2p22.2
Domains:	DUF185
MW:	49.2 kDa


[View online »](#)

Gene Summary:

This gene encodes an assembly factor protein which helps in the assembly and stabilization of Complex I, a large multi-subunit enzyme in the mitochondrial respiratory chain. Complex I is involved in several physiological activities in the cell, including metabolite transport and ATP synthesis. The encoded protein is a methyltransferase which methylates Arg85 of a subunit of Complex I in the early stages of its assembly. A pseudogene related to this gene is located on chromosome 8. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2016]