

Product datasheet for **RC202559A1V**

Human RPS25 (NM_001028) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human RPS25 (NM_001028) AAV Particle
Tag: Myc-DDK
Symbol: RPS25
Synonyms: S25
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC202559 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCCGCCTAAGGACGACAAGAAGAAGAAGGACGCTGAAAAGTCGGCCAAGAAAGACAAAGACCCAGTGA
ACAAATCCGGGGCAAGGCCAAAAAGAAGAAGTGGTCCAAGGCCAAAGTTCGGGACAAGCTCAATAACTT
AGTCTTGTGGTTCACAAAGCTACCTATGATAAACTCTGTAAGGAAGTTCCTCAACTATAAACTTATAACCCCA
GCTGTGGTCTCTGAGAGACTGAAGATTCGAGGCTCCCTGGCCAGGGCAGCCCTTCAGGAGCTCCTTAGTA
AAGGACTTATCAAAGTGGTTTCAAAGCACAGAGCTCAAGTAATTTACACCAGAAATACCAAGGGTGGAGA
TGCTCCAGCTGCTGGTGAAGATGCA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC202559 protein sequence
 Red=Cloning site Green=Tags(s)

MPPKDDKKKKDAGKSAKKDKDPVNKSGGKAKKKKWSKGVKVRDKLNNLVLFDKATYDKLCKEVPNYKLITP
AVYSERLKIIRGLARAALQELLSKGLIKLVSKHRAQVIYTRNTKGGDAPAAGEDA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_001028



[View online »](#)

ORF Size:	375 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001028.2</u>
RefSeq Size:	514 bp
RefSeq ORF:	378 bp
Locus ID:	6230
UniProt ID:	<u>P62851</u>
Cytogenetics:	11q23.3
MW:	13.7 kDa