

Product datasheet for **RC202561L4V**

DOK4 (NM_018110) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DOK4 (NM_018110) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DOK4
Synonyms:	IRS-5; IRS5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_018110
ORF Size:	978 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202561).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018110.2
RefSeq Size:	2750 bp
RefSeq ORF:	981 bp
Locus ID:	55715
UniProt ID:	Q8TEW6
Cytogenetics:	16q21
Domains:	PH, IRS
Protein Families:	Druggable Genome



[View online »](#)

MW: 37 kDa

Gene Summary: DOK proteins are enzymatically inert adaptor or scaffolding proteins. They provide a docking platform for the assembly of multimolecular signaling complexes. DOK4 functions in RET-mediated neurite outgrowth and plays a positive role in activation of the MAP kinase pathway (By similarity). Putative link with downstream effectors of RET in neuronal differentiation. May be involved in the regulation of the immune response induced by T-cells.[UniProtKB/Swiss-Prot Function]