

Product datasheet for **RC202567L3V**

EIF4A3 (NM_014740) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	EIF4A3 (NM_014740) Human Tagged ORF Clone Lentiviral Particle
Symbol:	EIF4A3
Synonyms:	DDX48; eIF-4A-III; eIF4A-III; eIF4AIII; Fal1; MUK34; NMP265; NUK34; RCPS
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_014740
ORF Size:	1233 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC202567).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_014740.2
RefSeq Size:	1734 bp
RefSeq ORF:	1236 bp
Locus ID:	9775
UniProt ID:	P38919
Cytogenetics:	17q25.3
Domains:	DEAD, helicase_C
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Spliceosome

MW: 46.9 kDa

Gene Summary: This gene encodes a member of the DEAD box protein family. DEAD box proteins, characterized by the conserved motif Asp-Glu-Ala-Asp (DEAD), are putative RNA helicases. They are implicated in a number of cellular processes involving alteration of RNA secondary structure, such as translation initiation, nuclear and mitochondrial splicing, and ribosome and spliceosome assembly. Based on their distribution patterns, some members of this family are believed to be involved in embryogenesis, spermatogenesis, and cellular growth and division. The protein encoded by this gene is a nuclear matrix protein. Its amino acid sequence is highly similar to the amino acid sequences of the translation initiation factors eIF4AI and eIF4AII, two other members of the DEAD box protein family. [provided by RefSeq, Jul 2008]